

INVACARE ACTION 3 COMFORT



NAME	INVACARE ACTION 3 COMFORT
PRODUCT CODE	CEG0041
SIZE	16", 18", 20"
FRAME	FOLDING POWDER COATED ALUMINIUM COLOUR: RED
BACKREST	RECLINING WITH GAS ASSISTED SHOCKS
FOOTREST	REMOVABLE ELEVATING
CASTOR	8X1" SOLID CASTOR
REAR RIM	24X13/8" ALUMINIUM SPOKE RIM WITH 24X13/8" GREY PU TYRE/REMOVABLE
UPHOLSTERY	GREY DARTEX COMFORT SEAT & BACK WITH HEAD REST
CARRY CAPACITY	LOAD CAPACITY – 125KG WEIGHT – 21KG
ARMS	HEIGHT ADJUSTABLE COMFORT ARMS
PACKING	108X92X42CM

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za