

BASIK PORTER

PRODUCT CODE: YK9070



The Basik Porter combines strength with thoughtful comfort features, making it the ideal choice for dependable, everyday mobility. Its 18" seat offers a perfect fit for most users, while height-adjustable armrests, swing-away footrests, and nylon upholstery ensure lasting comfort and support. Equipped with hub brakes for smooth, reliable stopping power, as well as quick-release 24" rear spoke wheels, lift wheels, and anti-tippers. It's built for both safety and ease of transport making it practical for caregivers. Delivering mobility, independence, and peace of mind in one sleek, durable package.



Hub Brakes



Anti Tipper



Double Cross Bar

Product Specifications




Width:
67cm




Seat Depth:
42cm



Seat Width:
18"



Rear Wheel:
24"



Front Wheel:
8"



Seat Height:
46cm



Length:
99cm



Height:
90cm



Back Height:
45cm



Max Weight Capacity:
135kg



Packaging DIM:
80 x 32 x 93cm



G / N Weight:
19.5 / 17kg

www.medop.co.za | @medopcc



Johannesburg:
011 827 5893
enquiries@medop.co.za
Medopark, 1 Manchester
Road, Wadeville

Cape Town:
021 534 8628
admin@medop.co.za
Viking Business Park, Corner
Park & Viking Road, Epping