

GLIDE RECLINER

PRODUCT CODE: YK9132



Designed for ultimate comfort and versatility, making it the perfect solution for users who need both mobility and extended seating support. With its reclining backrest ranging from 90° to 180° it offers personalized positioning for rest, recovery, or daily activities. The elevating leg rests and detachable desk-length armrests enhance comfort and accessibility, while the durable steel frame ensures long-lasting strength. Featuring 24" MAG wheels, anti-tippers, and hub brakes, the Glide Recliner prioritizes safety and ease of handling. Finished with nylon upholstery, this recliner wheelchair is ideal for both home and care facility use—delivering mobility with a touch of luxury.



Hub Brakes



Anti Tipper



MAG Wheels

Product Specifications




Width:
Seat width
+ 23cm




Seat Depth:
40cm



Seat Width:
16", 18", 20"



Rear Wheel:
24"



Front Wheel:
8"



Seat Height:
50cm



Length:
123cm



Height:
126cm



Back Height:
80cm



Max Weight Capacity:
120kg



Packaging DIM:
92 x 31 x 96cm



G / N Weight:
29.5 / 26.6kg



www.medop.co.za | [@medopcc](https://www.instagram.com/medopcc)



Johannesburg:
011 827 5893
enquiries@medop.co.za
Medopark, 1 Manchester
Road, Wadeville

Cape Town:
021 534 8628
admin@medop.co.za
Viking Business Park, Corner
Park & Viking Road, Epping