



CODE: YK9071S

The Glide SEMI-LITE is the latest folding frame model in the GLIDE range of wheelchair offerings. This attractive wheelchairs has numerous adjustments not normally found on an entry level model of wheelchair allowing the user to adapt the wheelchair for their own specific needs. The seat height adjustment is extremely useful for those users who require foot propulsion at an affordable cost and the lightweight spoke rims are removable for ease of transport. The Backrest TAB allows for positioning and this makes the Glide SEMI-LITE a favoured product for those users on a budget who require some level of configurability from their wheelchair.



- Frame: Steel
- Size: 16", 18", 20"
- Packing: 80x30x93cm N.W 18kg
- Upholstery: TAB
- Rear Wheel: 24 x 1 3/8" QR with solid PU tyre
- Castor: 7x2" Solid PU (compatible with 8x2")
- Adj: Adjustable Seat Height & Amputee Setting
- Load: 115kg

SIZE	16"	18"	20"
Overall Length (cm)	98	98	97
Overall Width (cm)	63	63	63.5
Overall Height (cm)	90	90	87
Overall Width Folded	33	33	31.3
Seat Width (cm)	41	46	51
Seat Depth (cm)	44	44	41.7
Back Height (cm)	42	42	45
Floor to seat Height	48	48	45
Weight Capacity (kg)	115	115	115
Front Caster Spec.	8"	8"	8"
Rear Wheel Spec.	24"	24"	24"
G./N. Weight (kg)	20.0/18.0	20.0/18.0	16.9/14.4
Carton Size (cm)	80x30x93	80x30x93	79x32x96



Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za