

GLIDE TRANSIT

PRODUCT CODE: YK9081T

NAPPI CODE:



Technical Specifications:

Seat Width	18"
Frame	Single cross brace, powder coating
Front wheel	8"*2" PU tires
Rear Wheel	Solid & QR 12" solid tyre
Armrest	Flip back slope armrest with PU armpad
Footrest	Swing in & away footrest
Seat and backrest	Nylon with sponge
Front fork	Plastic
Footplate	Plastic
Carton size	83x30x96cm
GW/NW	13.5kg / 11.5kg
Max. load capacity	120kg

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za