

PHOLA PRO POSITIONING WHEELCHAIR

Code: 12374-44



Specifications:

Width:	15", 17", 19"
Frame:	Alluminium
Arm Rests:	Detachable, Height Adjustable
Front Wheel:	8"x 32mm PU Tyre Caster
Rear Wheels:	24" Spoke wheel
Back Reclining Angle:	30 degrees
Seat Title Angle:	20 degrees
Footplate:	Angle adjustable footplate with one piece board
Leg Rests:	Swing in and out, detachable, elevating, calf pad
Accessories:	Anti Tip, Laterals, Pommel, Headrest
Packaging DIM:	73x22.5x86
Max weight Capacity:	150kg
∇.W:	12.3kg
G.W:	13.8kg

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za