

BASIK STEEL WHEELCHAIR

CODE: YK9041



Frame: Single cross brace, powder coating
Front Wheel: 8x1" Solid Castor
Rear Wheel: 24" spoke wheels, PU tire
Armrest: Flip back slope armrest, PVC armpad
Footrest: Swing away footrest, heel support
Seat and Backrest: Flame-retardant nylon with sponge
Front Fork: Plastic
Footplate: Plastic
Colour: Blue

Size 18"	SPECS MODEL YK9041													
		Overall length (cm) 103	Overall width (cm) 61	Overall height (cm) 91	Overall width folded (cm) 26.5	Seat width (cm) 46	Overall depth (cm) 40	Back height (cm) 43	Floor to seat height (cm) 50	Weight capacity (kg) 120	Front castor spec. 8"	Rear wheel spec. 24"	Carton size (cm) 82x28x93	NW/GW (kg) 21.5/19.0

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
 P.O. Box 51124, Raedene, 2124
 enquiries@medop.co.za | www.medop.co.za
 JHB Tel: 011 827 5893/4/5
 CT Tel: 021 534 8628 | admin@medop.co.za