

BASIK KIDDIES

The Basik Kiddies Wheelchair from Medop cc has been specifically designed with children in mind to meet a child's dimensions and requirements.

The wheelchair comes in a size appropriate 12" x 12" seat configuration for a child or for the young adult in a 14" x 16" and the arm rest heights and footrest have been designed with this in mind. The 6" solid front castor and 20" rear rim with solid PU tyre complete the wheelchair making these dimensions significantly easier for a child to propel as these sizes are anatomically more appropriate for optimal propulsion in these dimensions.

The wheelchair is completed with the new space saver flip back arm rest meaning an active child and busy parent can never lose the arm rest.

The Basik Kiddies boasts all the features one would expect to find on a higher end model designed to give the user even more from the wheelchair with

- Space Saver Flip Back Desk Arm Rests
- Wrap around reinforced back upholstery
- Fixed axle Spoke rims,
- Lightweight moulded nylon footplates
- Colour coded frame for appeal
- All packaged into a more manageable overall weight of 16kg.



*weight based on 12"x12" size

Specifications:

| | |
|--------------------|--|
| Frame: | Single cross brace, powder coating |
| Front wheel: | 6", PU tires |
| Rear wheel: | 20" Spoke wheel, PU tires |
| Armrest: | Flip back with desk length |
| Footrest: | Swing away footrest |
| Seat and backrest: | Nylon with sponge wrap around upholstery |
| Front Fork: | Aluminium |
| Footplate: | Plastic |

| SIZE | 12" | 14" |
|----------------------|-----------|-----------|
| Overall Length (cm) | 93 | 103 |
| Overall Width (cm) | 48 | 53 |
| Overall Height (cm) | 84 | 91 |
| Overall Width Folded | 27.5 | 26.5 |
| Seat Width (cm) | 31 | 35 |
| Seat Depth (cm) | 31 | 40 |
| Back Height (cm) | 38 | 41 |
| Floor to seat Height | 46 | 50 |
| Weight Capacity (kg) | 60 | 80 |
| Front Caster Spec. | 6" | 8" |
| Rear Wheel Spec. | 20" | 24" |
| G./N. Weight (kg) | 18.5/16.0 | 20.4/17.9 |
| Carton Size (cm) | 72x28x86 | 72x28x86 |