



Code: YK9080H

The Lightweight Glide Lite has awed the market with the ability to adjust for a comfortable seat depth for a standard user all the way up to a 50cm seat depth for the unusually tall user. Packaged in a lightweight, user friendly aluminium folding frame the Glide Lite allows for arm rest height adjustment to compliment the taller user, angle adjustable footplates on a reinforced triple-cross-bar and a TAB for positioning. There is not much you can't achieve in seating, functionality and comfort with the Glide Lite



- Frame: Aluminium
- Sizes: 15", 16", 18", 20"
- Castor: 8x2" Solid PU
- Rear Rim: 24x1 3/8" QR rim with Solid PU Tyre
- Arm: Flip Back & Height Adj Sport Arm
- Footrest: Swing Away and removable with height adjustment
- Footplate: Angle Adjustable
- Upholstery: Tension Adjustable Nylon Upholstery
- Adjustment: Adjustable COG & Seat Angle
- Back: Height Adjustable with Push Handle
- Seat Depth: 42cm to 51cm
- Packing: 84x29x78cm
- N.W: 15.5kg
- Load Capacity: 115kg

SIZE	15"	16"	18"	20"
Overall Length (cm)	101	101	101	101
Overall Width (cm)	59	61.5	65.5	70.5
Overall Height (cm)	92	92	92	92
Overall Width Folded	32	32	32	32
Seat Width (cm)	38	41	46	51
Seat Depth (cm)	41-51	41-51	41-51	41-51
Back Height (cm)	48.5	48.5	48.5	48.5
Floor to seat Height	50	50	50	50
Weight Capacity (kg)	115	115	115	115
Front Caster Spec.	8"	8"	8"	8"
Rear Wheel Spec.	24"	24"	24"	24"
G./N. Weight (kg)	16.8/14.8	17.3/15.3	18.3/16.3	19.3/17.3
Carton Size (cm)	84x29x78	84x29x78	84x29x78	84x29x83



Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za