



Code: YKCOMFORT

Users requiring almost all the benefits of an Active wheelchair in a folding frame format need look no further than the Glide PERFORMANCE Wheelchair from Medop cc. This wheelchair has just about all the adjustability an active user will require in a lightweight folding wheelchair configuration. The TAB, angle adjustable footplate, Adjustable seat depth, 6" performance castor and bent over mud-guards all packaged in a sleek, aesthetically pleasing design make the Glide PERFORMANCE hard to match and fits into a wheelchair category all of its own.



- Frame: Aluminium
- Sizes: 15", 16", 17", 18'
- Castor: 6x1.25" Solid
- Rear Rim: 24x1 3/8" rim with Pneumatic Tyre
- Arm: Removable Half-Arm with foam
- Footrest: Swing Away and removable with height adjustment
- Footplate: Angle Adjustable
- Upholstery: Tension Adjustable Nylon Upholstery
- Adjustment: Adjustable COG & Seat Angle
- Back: Height Adjustable with Push Handle 81-97cm
- Seat Depth: 42cm
- Packing: 80x32x86cm N.W: 13.5kg
- Load Capacity: 115kg

SIZE	15"	16"	17"	18"
Overall Length (cm)	93	93	93	93
Overall Width (cm)	55.5	58	60.5	63
Overall Height (cm)	81-97	81-97	81-97	81-97
Overall Width Folded	30	30	30	30
Seat Width (cm)	38	41	43	46
Seat Depth (cm)	42	42	42	42
Back Height (cm)	36-37	36-37	36-37	36-37
Floor to seat Height	46-52	46-52	46-52	46-52
Weight Capacity (kg)	115	115	115	115
Front Caster Spec.	6"	6"	6"	6"
Rear Wheel Spec.	24"	24"	24"	24"
G./N. Weight (kg)	14.5/12.5	15.0/13.0	15.5/13.5	16.0/14.0
Carton Size (cm)	80x32x86	80x32x86	80x32x86	80x32x86



Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za