



Code: YKACTIVE

Active wheelchair users will be impressed with the ease of manoeuvrability and propulsion provided by the Glide SPRY Active Wheelchair. There is the option for an anti-tipper for the new wheelchair user or Off-Road rims for the more rural user. The Glide SPRY also offers user the unique option of a front fold-up footplate to allow users improved toilet access. Medical professionals will be impressed with the numerous seating and positioning features available on the Glide SPRY and the wheelchair was specifically designed to allow users increased freedom within the Compensation Fund's budget allowance as well as various other funders Above Threshold External Medical Appliance benefits.



**Features:**

- Aluminium Powder Coated Frame
- 5x1.5" front castors
- QR Light Weight pneumatic rear spoke wheels
- Adjustable COG
- Single piece footplate. Height Adjustable
- Seat and Backrest – Black Nylon
- Seat Depth: 16" , 18"
- TAB back upholstery
- Height Adjustable push handle
- Flip down & Angle adjustable backrest
- Undermount brake system
- Colours: Black & Red / Red & Blue
- Option: Anti-Tipper

| SIZE                 | 16"       | 18"       |
|----------------------|-----------|-----------|
| Overall Length (cm)  | 79        | 79        |
| Overall Width (cm)   | 64        | 69        |
| Overall Height (cm)  | 82-86     | 82-86     |
| Overall Width Folded | /         | /         |
| Seat Width (cm)      | 40.5      | 45.5      |
| Seat Depth (cm)      | 41        | 41        |
| Back Height (cm)     | 33        | 33        |
| Floor to seat Height | 47-51     | 47-51     |
| Weight Capacity (kg) | 120       | 120       |
| Front Caster Spec.   | 4"/5"     | 4"/5"     |
| Rear Wheel Spec.     | 24"       | 24"       |
| G./N. Weight (kg)    | 14.0/12.0 | 14.0/12.0 |
| Carton Size (cm)     | 68x62x60  | 68x62x60  |



Medopark, 1 Manchester Road, Wadeville,  
 P.O. Box 51124, Raedene, 2124  
 enquiries@medop.co.za | www.medop.co.za  
 JHB Tel: 011 827 5893/4/5  
 CT Tel: 021 534 8628 | admin@medop.co.za