



Code: YK9063

The Glide SWIFT has proven to be the most popular wheelchair model for post CVA rehabilitation discharge. With the reinforced double crossbars, height adjustable arms, TAB, Quick Release lightweight rims all mounted on a lightweight frame for transport, the wheelchair is manageable for elderly couples to handle. The robust and reinforced frame have also made this wheelchair a favourite with users looking for a durable primary or back-up wheelchair on a reasonable budget. Therapists will be enamoured with the user friendliness of the Glide SWIFT in confidently allowing the therapist to prescribe their patients with a clinically appropriate wheelchair for discharge.



TAB UPHOLSTERY



ADJUSTABLE T-ARM



FRAME	Double cross braces, liquid painting
FRONT WHEEL	8"x1.25", PU tires
REAR WHEEL	24", QR spoke rim with PU tyre
ARMREST	Flip back and Height Adjustable T-arm
FOOTREST	Swing in & away footrest
SEAT AND BACKREST	Nylon with sponge
FRONT FORK	Aluminium
FOOTPLATE	Plastic

SIZE	16"	18"	20"
Overall Length (cm)	103	103	103
Overall Width (cm)	60	65	70
Overall Height (cm)	92	92	92
Overall Width Folded	30	30	30
Seat Width (cm)	43	46	51
Seat Depth (cm)	41	41	41
Back Height (cm)	43	43	43
Floor to seat Height	49	49	49
Weight Capacity (kg)	135	135	135
Front Caster Spec.	8"	8"	8"
Rear Wheel Spec.	24"	24"	24"
G./N. Weight (kg)	15.8/13.8	16.3/14.3	16.8/14.8
Carton Size (cm)	82x31x93	82x31x93	82x31x93



Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za