

# JADE PRO

PRODUCT CODE: WCX4



## Heavy Duty Light Weight Wheelchair

- Height adjustable back
- Angle adjustable back
- Push handle height adjustable
- Tension adjustable backrest upholstery
- Flip back and detachable armrest with height / length adjustable arm pad

### Specifications:

Seat Width :	22" / 24" / 26" / 28" / 30"
Seat Depth :	Adjustable. Starting from 44cm to 50cm
Load Capacity:	205kg
	Multiple COG positions
Front wheel:	200 x 50 mm solid PU caster
Back wheel:	QR 24x 1 3/8" spoke HD wheel
N.W:	28kg
G.W:	34kg
Ctn Size:	83 x 39 x 105 cm

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za