

LULA LITE

Code: CA9881LQ

Active users will not only thoroughly enjoy the **Lula Lite** sleek look but once a user sits in the **Lula Lite** they will immediately notice the impressive performance of this wheelchair. With Carbon Fibre Mud Guards and a Carbon Fibre Axle Bar this super lightweight wheelchair provides the Active user increased independence without compromising on performance. The **Lula Lite** is designed with the Active user in mind and multiple options such as arm rest, anti-tipper and off-road rims make this a wheelchair of preference for any Active user.



Alu Castor



Brakes & Spokes



Mudguard



Footplate

Seat Width:	15", 16", 17", 18"
Seat Depth:	17"
Frame:	80° Reinforced Aluminium with 8° seat and Carbon Fibre Axle Bar
Seat:	Active Light Tension Adjustable Seat
Back:	Active Light Tension Adjustable Back
Back Height:	32cm-41cm
Back Angle:	87° - 93°
Side Panel:	Carbon Fibre Mud Guards
Rims:	24" Colour Coded Performance Spoke Rims
Tyres:	24x1" Pneumatic Tyre
Castor:	Aluminium 5x1 ¼" solid Castor
Brake:	NEW Active Lightweight Brake
Footplate:	Angle Adjustable
Adjustment:	Centre of Gravity Adjustment
Weight:	GW11.5kg/NW9.8kg
Transport weight:	7.6kg
Load Capacity:	115kg
Options:	Anti-Slip Push-Rim Covers, Bolt-On Arm Rest, Anti-Tip, Off-Road Rims, Seat Bag

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za