

SATURN HERCULES^{NG}

The New Generation (NG) Saturn Hercules from Medop cc has all the great qualities you have come to expect from the Saturn Range of durable, economic steel frame folding wheelchairs including the 8x2" solid castor and 24x1.75 PU rear tyre.

The New Generation of Saturn wheelchairs have been designed to give the user even more from the wheelchair with

- Space Saver Flip Back Arm Rests making it easier to go through a doorway,
- Wrap around reinforced back upholstery,
- Quick Release Rear MAG rims,
- Lightweight moulded nylon footplates

All packaged into a more manageable overall weight of 18kg.

The incredible Saturn Hercules NG is now available from 14" all the way through to a Bariatric 24" model with a load capacity from 120kg up to 150kg and the **double reinforced crossbars** and frames on our 20" sizes upwards.



* supplied with MAG rims
*weight based on 16"x16" size

SIZE	14"	16"	18"	20"	22"	24"
Overall Length (cm)	103	103	103	103	103	103
Overall Width (cm)	53	58	63	68	73	78
Overall Height (cm)	91	91	91	91	91	93
Overall Width Folded	26.5	26.5	26.5	26.5	26.5	26.5
Seat Width (cm)	35	40	45	51	55	60
Seat Depth (cm)	40	40	40	40	40	40
Back Height (cm)	41	41	41	41	41	46
Floor to seat Height	50	50	50	50	50	50
Weight Capacity (kg)	80	120	120	135	150	150
Front Caster Spec.	8"	8"	8"	8"	8"	8"
Rear Wheel Spec.	24"	24"	24"	24"	24"	24"
G./N. Weight (kg)	20.4/17.9	20.6/18.1	21.0/18.5	22.5/20.0	23.0/20.5	23.5/21.0
Carton Size (cm)	72x28x86	81x27x94	81x27x94	81x27x94	82x27x93	82x28x96

Specifications:

Frame:	Single cross brace, powder coating
Front wheel:	8x2" solid castor
Rear wheel:	24" Quick Release Mag rim with 24x1.75" PU tyre
Armrest:	Flip back with desk length
Footrest:	Swing away footrest
Seat and backrest:	Nylon with sponge
Front Fork:	Steel
Footplate:	Plastic



Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za