


Aquatec | Shower stools, shower chairs



Invacare I-Fit 9781E



Invacare I-Fit 9781E

- Height Adjustable from 38 cm to 55 cm
- Removable backrest
- Ergonomic shaped seat with backrest for extra comfort
- Easy to clean
- Light weight and easy to transport



Technical data

Order numbers	Aquatec Pico
Shower stool	9.20.001
Shower stool with armrests	9.21.001
Shower chair with backrest	9.23.001
Shower chair complete	9.22.001
Backrest	16469
Armrests (pair)	16468

Dimensions	
Seat width	43 cm (A)
Seat depth	42 cm (B)
Seat height	42,5 - 57,5 cm (C)
Total width	57,5 cm (D)
Total depth	52 cm (E)
Total height	82,5 - 97,5 cm (F)
Height to armrests	60,5 - 75,5 cm (G)
Width between the armrests	43 cm (H)
Width, hygiene recess	18 cm
Depth, hygiene recess	14,5 cm
Graduation, height adjustment	2,5 cm

More info	
Weight, shower stool	3,2 kg
Weight, shower stool with armrests	4,1 kg
Weight, shower chair with backrest	4,0 kg
Weight, shower chair complete	4,8 kg
Load capacity	160 kg!
Colour	White

Technical data

Order numbers	Invacare I-Fit 9781E
Shower stool	
Shower chair	9781E

Dimensions	
Seat width	34
Seat depth	37
Seat height	38-55
Total width	50
Total depth	43
Total height	76-93

More info	
Weight	4,5 kg
Load capacity	180 kg
Colour	white

AQUATEC



Yes, you can!

Medopark, 1 Manchester Road, Wadeville, P.O. Box 51124, Raedene, 2124

medopcc@global.co.za, www.medop.co.za

Tel: 011 827 5893/4/5

Fax: 011 827 6286

pg. 35

