

DOLOMITE | *STEP UP*



The rehab optimizer

DOLOMITE STEP UP is a new modern platform support. Made to give rehabilitation to several patients every day. Adjust to optimal ergonomic position by a few simple handgrips – quick and easy!

DOLOMITE STEP UP can be manually adjusted, hydraulically or with a motor. Assembly of accessories is tool free. Are you ready to step up?





Made for training

The base frame is designed to enhance premium gait area. The raised frame in front allows longer step. The width in side the frame generous.

Ergonomic adjustment of arm support

Arm pads are in soft material for pressure relief. Securely assembled with microfilm. Quick and simple to remove for cleaning or replacement.

Transfer and access

Well thought out for transfer situations. From bed, the lowered rear-end helps reversing under the bed. When

transferring from wheel chairs benefit of the generous frame size and, the width adjustment.

Features and Options



By releasing the central knob, the both arm pads are simultaneously adjusted in width. The position of the arm supports remains in an angle of 25 degree from the shoulder points - the ideal position for keeping upper strength.

Accessories and options

Quick and easy assembly in track or prepared recess.



Practical basket

Keeping belongings within easy to reach range.
Art.no 1531390



Tray

for small things.
Art.no 1531391



Elbow stop

Securing elbows from falling out of position.
Art.no 1531448



Swivel caster with directional brakes

Helps focusing on walking instead of direction.
Art.no 1531408



Brake system

For the user capable of operating brake function.

Hand-brake kit

Art.no 1537656

One-hand brake

Art. no Right 1539513
Art. no Left 1531450



Oxygen tank holder

Secure mounting of oxygen tube on the walker.
Art.no 1531406



IV holder

Separate IV-stand makes life complicated. Not anymore.
Art.no 1531407



Stand assist inkl. slings

For verticalisation. Fits hydraulic and manual version.

Art.no 1531449 (standing plate + knee support + sling)

Standing plate

Art.no. 1535330



Ball grips

when having difficulties to close the hand.

Art.no 1539515

Others

Small caster 100 mm diameter.
Art. no 1536047

Models



Manual
Pull out the pin move to required position release the pin.



Pull up Hydraulic
Pull up lever to adjust to required height.



Electric
Using the hand control to adjust to required height. The motor version gives active assistance.

Technical data



Step Up Manual	580/670/760 mm	700/790/880 mm	870 - 1320 mm	340 - 470 mm	620 - 750 mm
Step Up Hydraulic	580/670/760 mm	700/790/880 mm	870 - 1320 mm	340 - 470 mm	620 - 750 mm
Step Up Electric	580/670/760 mm	700/790/880 mm	870 - 1320 mm	340 - 470 mm	620 - 750 mm



Step Up Manual	810 mm	16 kg	150 kg	125 mm	215 mm
Step Up Hydraulic	810 mm	18 kg	150 kg	125 mm	215 mm
Step Up Electric	810 mm	23 kg	150 kg	125 mm	215 mm

Colour



RAL 1035

Article code

Step Up Manual 1521677
Step Up Hydraulic 1521678
Step Up Electric 1521679

Invacare Dolomite AB
Växjövägen 303 75
Box 200, 343 75 Diö
Sweden
Tel. +46 (0)476 535 00
Fax +46 (0)476 535 99
www.dolomite.biz

© 2011 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

StepUp - EU - 08/2011

