

INVACARE ADEO TOILET FRAME



*Maris Toilet Seat Raiser for Illustration Purposes

NAME	Invacare Adeo Toilet Frame
Frame	Aluminium
Base	550 x 440 (mm) on adjustable feet for stability
Height	Adjustable 650 to 750 (mm)
Carry Capacity	Weight – 100kg Gross Weight – 2.5kg
Packing	670 x 440 x 100(mm)



Medopark, 1 Manchester Road, Wadeville, P.O. Box 51124, Raedene, 2124

medopcc@global.co.za, www.medop.co.za

Tel: 011 827 5893/4/5

Fax: 011 827 6286