

ADJUSTABLE BATH BOARD



2) Slatted Base

NAME	Adjustable Bath Board
Seat	Adjustable one piece moulded plastic
Packing	750 x 340 x 140(mm)
Seat Size	700 x 310 (mm)
Weight	2.8kg Gross Weight 3.3kg
Load Capacity	100kg
Options	1) With TRANSFER Disc 2) Slatted Base



Medopark, 1 Manchester Road, Wadeville, P.O. Box 51124, Raedene, 2124

medopcc@global.co.za, www.medop.co.za

Tel: 011 827 5893/4/5

Fax: 011 827 6286