

# ALUMINIUM CABLE BRAKE WALKER



## Features:

- Lightweight Aluminium
- Height adjustable handle
- Light Squeeze hand brakes with static lock function
- 8" solid castors
- Padded Vinyl seat
- Large storage basket



Width inside handgrips:	48cm
Overall width:	61cm
Overall depth:	54cm
Seat area:	35x32x3.5cm
Seat height:	49cm
Adjustable handle height:	74-85.5cm
Front wheel:	6 inch (PVC)
Rear wheel:	8 inch (PVC)
Packing dimension:	76x23x61cm
Max user weight:	136kg
Weight:	6.7kg

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za