

BARIATRIC H/DUTY ELECTRIC BED



Function:	Electric Hi-Lo 3 Function
Castors:	5" PA Castors (foot w/brake)
Frame:	Reinforced Steel – Powder Coated
Outer Dimensions:	2225x1070mm
Surface Dimension:	2020x880mm
Adjustable Height:	460-960mm
Bed Head Adj:	To 70 degrees G.W: 136kg
Packing:	Foot: 117x23x112cm (53kg) Head: 112x23x110cm (42kg) H&F Board: 112x19x70 (34kg)
Upper Side Rail	74x7x45cm (7kg)
Includes:	Trunk Side Rail Set, Emergency Manual Crank
Load Capacity:	247kg

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za