

CLASSIC VISCO CUSHION

When simplicity and effectiveness work together

The Invacare® Classic Visco cushion is entirely composed of viscoelastic shape memory foam with a density about 85 kg/m³. There are two firmnesses: firm and very firm, depend on the sizes. With its height of 75mm, this flat cushion Classic Visco is designed to aid prevention of skin breakdown with medium to high risk levels. It is delivered with two covers.



- 100% viscoelastic foam
- Flat cushion
- Two-way elastic, breathable, waterproof fabric cover
- 4 sizes available

Features



A waterproof outer cover, for a tight cushion



Handle

Located in a corner, this handle is useful and discreet!



Two firmnesses, 4 sizes

Firm or very firm, depending on the size.

Technical data



Classic Visco

Width :
375 mm
405 mm
430 mm
455 mm



Depth :
400 mm
430 mm



Hight :
75 mm



From
0,950 kg



Firm : Max. 50
kg

Very firm : max
100 kg



Outer cover
cleaning : 90°

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za