

# CAPE WALKER

The Cape Walker was designed to ease the learning and relearning of walking. The locking system on the rear wheels will prevent the walker from moving backwards and will encourage the user to move forward. It also has a directional front wheel setting.

The system is fully reversible so that your walker can be used in a forward direction as well and is foldable for transport and storage.



## Technical Information:

- Rounded frame/caster Joint
- Foldable
- Hip Aligners
- Comfort Support Sling seat
- Anti-reverse Lock
- Directional Wheel lock

Code	Colour	Height adjustment
11261	Gold	37.5 - 47.5cm
11262	Bronze	47.5 - 65cm
11263	Blue	60 - 80cm
11264	Silver	72.5 - 92.5cm

## How to pick your walker

It suffices to pick a Cape Walker where the height of the handles is equal to the height of the user. Every walker is adjustable in height for a perfect personal solution.