

CLOUD CAR SEAT



The Cloud Car Seat is meticulously crafted to redefine every car journey for children from 3 to 12 years old with unique requirements. The Cloud car seat offers multiple seating modes which can meet the child's changing needs during growth. It can also integrate with different optional accessories, allowing for personalized customization that guarantees utmost safety and supports the user in maintaining an optimal and cozy posture throughout their road adventures. The Cloud car seat embodies safety, portability, and easy assembly, offering both the user and caregiver a delightful and hassle-free car voyage.



Adjustable Headrest



Adjustable Side Wings




Convenient to Carry

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

CLOUD CAR SEAT



Adjustable Seating Modes

The headrest's height is able to be adjusted as the user grows. The height and width of the side wings are adjustable, providing a comfortable and stable seating experience for users of different body shapes.

Utmost Safety

The Cloud car seat comes equipped with a five-point harness system that secures the child firmly in place while also accommodating safety belts and any additional positioning harness they may require. The harness is adjustable to ensure a snug fit without causing discomfort.

Personalized Customization

We offer a variety of optional accessories, such as side wings and leg rests, to meet the unique needs of different users, enhancing both their comfort and safety during travel.

Quick Installation

The Cloud car seat features the ISOFIX connect system, which makes it fit into a wide range of vehicles, and enables easy and rapid installation and removal. Also, the Cloud is lightweight and foldable, allowing for storage in the trunk.

Accessories

All seat cushions and accessories are crafted from breathable, antibacterial materials, allowing for effortless removal and cleaning.



Calf Stabilizer

The calf stabilizer is designed for the user with calf muscle weakness, instability or require extra support while sitting. It can help keep the user's calves in a fixed position, effectively preventing crossing or "scissoring".



Side Wings

The side wings offer upper body support, ensuring proper posture and enhancing both comfort and safety during car journeys.



Item Dimensions	1
Headrest Width (A)	31 cm
Overall Width (B)	40 cm
Overall Length (C)	36 cm
Overall Height (D)	70cm
User Age	3 - 12 years
Frame Weight	6 kg
Weight Capacity	36 kg