

CONTOURA CUSHION

PRODUCT CODE: CON



This contoured memory foam cushion is designed to provide enhanced comfort and pressure redistribution for permanent users. The high-quality memory foam moulds to the body, helping to support proper posture and improve seating stability. Its ergonomic contoured design promotes pressure relief and increased comfort during extended periods of sitting. The cushion includes a protective removable cover for easy cleaning and durability. Available in various widths to suit different seating requirements.



Removeable
washable cover




Memory Foam




Contoured for
pressure relief


Product Specifications



Depth:
16"



Thickness:
8.5cm



Width:
16", 18", 20", 22"



Cover Wash
Temp:
30 degrees



Max Weight Capacity: Product Weight:
110 / 127kg 42x30x61cm



www.medop.co.za | @medopcc



Johannesburg:
011 827 5893
enquiries@medop.co.za
Medopark, 1 Manchester
Road, Wadeville

Cape Town:
021 534 8628
admin@medop.co.za
Viking Business Park, Corner
Park & Viking Road, Epping