


# CUDDLE POSITIONING SYSTEM

Code: BSK002 / BSK063

**BSK002**  
Includes:  
Base Mattress M  
Horseshoe Pillow  
Lateral Support Pillow  
Leg support Pillow  
Calf stabiliser  
Chest strap  
Storage bag


**BSK063**  
Includes:  
Base Mattress L  
Body Pillow  
Foot support pillow  
Cylindrical pillow  
Chest Strap  
Neck support Pillow




The cuddle positioning system is designed specifically for individuals who require special assistance during lying or sitting posture training, providing them with a comfortable resting and training environment. The mattress is available in two sizes to accommodate children during different growth stages

## Base Mattress

The base mattress incorporates advance technology to alleviate pressure points and preventing bedsores. The mattress is crafted from hypoallergenic materials. The cover of mattress features two different breathable fabrics that helps regulate body temperature and moisture and is able to be removed and washed by machine, ensuring a clean and hygienic sleep environment.



Item No	A	B
BSBM02	160cm	65cm
BSBM03	190cm	65cm



Item No	A	B	C
BSWP01	48.5cm	24cm	24cm


## Wedge pillow

The wedge pillow is designed with a unique triangular shape that offers multiple benefits for various sleeping positions and health conditions.

Its gradual incline provides gentle elevation, which can provide support to the lumbar spine by placing it behind the back or improve leg blood circulation and alleviate lower back pressure by placing it behind the legs.

## Horseshoe pillow


The horseshoe pillow provides exceptional neck and head support. The special fabric on the headrest of the pillow provides breathability while preventing sweat and other fluids from sleeping into the pillow.



Item No	A	B
BSHP01	36cm	50cm

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za




Item No	A	B	C
BSBP01	62cm	60cm	20cm

### Body Pillow

The u shape design of the body pillow allows for versatile positioning.


The user can wrap their arms and legs around it to relieve pressure points and support their back.



Item No	A	B	C	D
BSNS01	45cm	21.5cm	10cm	22cm

### Neck support pillow

The contoured design of the neck support pillow can help fix the users head position and promotes proper alignment of the user's cervical spine, helping to alleviate tension and stiffness.




Item No	A	B
BSCP01	56cm	13cm

### Cylindrical Pillow


The cylindrical pillow can be used as a bolster pillow to support neck, lower back or legs, providing customized comfort for various sleeping positions.

### Lateral support pillow

The lateral support pillow is designed to help users to stabilize their body position while lying down, maintaining a correct and comfortable posture, relieving pressure on the spine and muscles and preventing users from rolling or falling.



Item No	A	B	C
BLS01	29cm	8cm	17cm



Item No	A	B	C
BSLP01	23cm	12.5cm	12.5cm


### Legs support pillow

The leg support pillow can elevate the users legs, relieving pressure on joints and muscles and promoting better blood flow.


### Calf stabilizer

The calf stabilizer is specifically designed to help individuals who experience calf muscles weakness, instability or require extra support while lying.

It can help keep the users calves in a fixed position, effectively preventing crossing or scissoring.



Item No	A	B	C	D
BSCS01	33cm	15cm	14cm	12cm



Item No	A	B	C	D
BSFS01	20cm	17.5cm	11cm	14cm

### Foot support pillow

The foot support pillow helps maintain proper alignment of the feet ankles and lower limbs while lying down. It can help relieve pressure on heels, ankles and soles.

