

DELTA THREE WHEEL WALKER



Features:

- Steel powder coated frame
- 8" PVC wheels provide more comfort
- Includes shopping basket and bag
- Push-down braking system
- Easy folding system for storage



Overall width:	64cm
Overall depth:	59cm
Tray height:	70cm
Adjustable handle height:	82-94cm
Front wheel:	8 inch (PVC)
Rear wheel:	8 inch (PVC)
Packing dimension:	66x31x71cm
Max user weight:	136kg
Weight:	7.6kg

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za