

DUDUZA ELECTRIC HOME CARE BED



Product Code: 9016-2

Economy Electric Homecare Bed

- Steel frame, surface with coating
- Bed deck with 4 sections
- Anti-sliding steel surface
- Density wooden headboard/footboard
- Aluminium side rails, coloured like beech wood
- 4" casters with brake
- Moteck 3 motors and 5 functions hand controller
- Bed Surface Size: 2000mm x 900mm
- Overall Length: 2120mm
- Overall Width: 1030mm
- Bed Surface Height: 250-750mm
- Max. Inclination of Back Section: 70°
- Max. Inclination of Leg Section: 25°
- Max. Inclination of Overall Bed: 12°
- Input Voltage: 220V50HZ
- Weight Capacity: 180kg
- N.W.: 106 kg, G.W. 118 kg
- Packing: 1/3ctns (ctn size: 930 x 390 x 1240mm)
- **Option: Backup battery, lifting pole, IV pole.**

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za