

ELECTRIC HI-LO HOSPITAL BED

Code: ALK06-B02P



4 sections

Specifications:

Width:	98cm
Length:	213cm
Load Capacity:	150kg
Batteries:	Back up unit standard
Actuators:	3 Electric
Hi-Lo:	45cm – 70cm
Back Rise:	80 degrees
Foot Lift:	40 degrees at the knee
Material:	Powder Coated Mild Steel
Head & Foot:	Removable ABS Head & Foot Board
Castors:	5" Dust Free rubber castor with brake on 2 castors
STD Accessories:	Cotsides, 12cm Taffeta covered Mattress, I.V Pole

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za