

ENGULM POWER WHEELCHAIR

PRODUCT CODE: ENGULMPOWER18G




The Engulm Electric Wheelchair, an impressive rear wheel drive wheelchair made from a powder coated steel frame


Designed with the following standard features: tilt, recline, headrest, electromagnetic brakes, lights, height adjustable armrests, detachable and width adjustable footrests with heel support



Lights



14" Pneumatic Rear Wheel



Headrest

Product Specifications




Width:
62.5cm




Seat Depth:
33 - 50cm



Seat Width:
45cm



Rear Wheel:
14" pneumatic



Front Wheel:
9" pneumatic / PU Tyre



Back Height:
49.5cm



Length:
114.3cm



Height:
99.5cm



Packaging DIM:
102 x 77 x 97cm



Max Weight Capacity:
135kg



Battery:
24V 50AH



Gross / Net Weight:
110kg / 98kg

Medop

www.medop.co.za | @medopcc



Johannesburg:
011 827 5893
enquiries@medop.co.za
Medopark, 1 Manchester Road,
Wadeville

Cape Town:
021 534 8628
admin@medop.co.za
Viking Business Park, Corner Park &
Viking Road, Epping