

GEZA COMMUNE RANGE

AQUA Taxi Shower Commode

AQUA18 – 5



Product Code: CA6204L

- Aluminium powder coated frame
- PU seat
- PU armrest & backrest
- Flip up armrest
- Swing away footrest
- 5" solid castors
- Foldable frame
- Strong cross bar foldable frame

Width 55cm	Seat Width 46cm	Seat Depth 40cm	Seat Height 51cm	Length 89cm	Back Height 45-48cm	Height 95-98cm	Max Load 120Kg	NW/GW KGS 16/18	Packaging 70x28x74cm	PCS/CTN 1

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za