



PATIENT LIFT

Code: Y401



FRAME	Steel, powder coating
CHASSIS LENGTH	1150mm
CHASSIS WIDTH	580mm - 1080mm
LIFTING HEIGHT	870mm - 1845mm
HANDLE HEIGHT	930mm
CASTOR HEIGHT	95mm
MOTOR	DC 24V 8000N
BATTERY	DC 24V, 5AH, Approx. 60 lifts per full charge
CHARGER	24V-1A
CHARGING TIME	6 hours per time
FRONT CASTOR	3" solid tires
REAR CASTOR	5" solid tires with locking
PUSH HANDLE	Push handle with rounded corners
MAX. USER WEIGHT	150KGS (23.62 stone)
SLING	Purpose designed comfort sling
PACKAGE	1240×630×360mm
G./ N. W.	50 / 46.5KGS



Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za