



VIRGO-D SCOOTER

Code: VIRGO-D

Meet your ideal Travel Companion in the Virgo-D. Whether using the Virgo-D for a sunset ride on the promenade or to meet up with your friends in your estate, the Virgo-D will get you there in style and the quick disassembly makes the Virgo-D easy to take along on your road trip adventures.



Frame:	Steel and plastic cover
Motor:	MTM 270W
Max Speed:	8km/h
Battery:	Lead acid 12AH
Travel Mileage:	10km
Controller:	Dynamic DR-50
Slope:	9 degree
Front Wheel:	8" black PU tyres
Rear Wheel:	8" black PU tyres
Armrest:	Flipback armrest
Seat:	PU / ABS + soft sponge, forward and backward adjustable
Special:	Easy to disassemble for storage and transfer
Max Capacity:	136kgs



Medop cc


Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za



VIRGO-D SCOOTER



Dimensions (cm/in)




Unfold



Weight(kg/lb)

N.W/Batt.	Battery/Ea	G.W/Batt.	Capacity
45/20.5	1.7/3.7	49/22.5	136/62



N.W. ≤ 45.0 kg
 G.W. ≤ 49.0 kg
 PACKING SIZE 1070x560x600 mm
 LOADING DETAILS (FOR REFERENCE ONLY)



Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
 P.O. Box 51124, Raedene, 2124
 enquiries@medop.co.za | www.medop.co.za
 JHB Tel: 011 827 5893/4/5
 CT Tel: 021 534 8628 | admin@medop.co.za