

GUTTER CRUTCH

Specifications:

- Available in Adult and Kiddies
- Adult adjustable from 95cm to 118cm
- Kiddies adjustable from 75cm to 98cm
- Arm pad length 23cm – padded with forearm strap
- 19mm reinforced ferrel
- Grey Colour



Maximum
User Wt:



Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za