



FIXI

Aluminium rollator for children

Equipment:

- dual action brake
- foldable
- airless tyres
- particular light powder coated aluminium frame
- approved up to 100 kg body weight

ref. no. 186.00.10

Accessories:

- seat [ref. no. 286.65.25](#)
- back pad [ref. no. 280.10.48](#)

Only 6,2 kg!



ref. no. 186.00.10



Rollators

9



ref. no. 186.00.20



FOX

Aluminium rollator for children/youth/ladies

Equipment:

- dual action brake
- foldable
- airless tyres
- particular light powder coated aluminium frame
- approved up to 100 kg body weight

ref. no. 186.00.20

Accessories:

- seat (black) [ref. no. 286.65.25](#)
- back pad [ref. no. 280.10.48](#)

Technical data:

Type:	FIXI	FOX
depth:	68 cm	68 cm
width:	58 cm	58 cm
height:	60,5 cm - 69,5 cm	69,5 cm - 77,5 cm
seat height:	46,5 cm	46,5 cm
seat width:	38,5 cm	38,5 cm

Type:	FIXI	FOX
weight:	6,2 kg	6,3 kg
adjustments in height:	5 x	5 x
max. body weight:	100 kg	100 kg
ref. no.:	186.00.10	186.00.20

Technical changes and changes in colour reserved.