


HINGE GRAB RAIL


PRODUCT CODE: H330/1



The Basica Hinge Grab Rail is a durable foldable support rail designed to provide extra stability and confidence in the bathroom. Manufactured from double steel tubing with a corrosion-resistant plastic coating, it offers reliable support for transfers and toilet use. The wall-mounted folding design saves space when not in use, making it ideal for both home and clinical environments. With a maximum user weight of 105 kg, the Hinge Grab Rail is a practical and dependable solution for users with limited mobility.



Foldable



Corrosion
Resistant



Wall Mounted

Product Specifications



Max User Capacity:
105kg



Frame Colour:
White



Product Weight:
2.2kg



Packaging DIM:
10x67x20.5cm



www.medop.co.za | @medopcc



Johannesburg:
011 827 5893
enquiries@medop.co.za
Medopark, 1 Manchester
Road, Wadeville

Cape Town:
021 534 8628
admin@medop.co.za
Viking Business Park, Corner
Park & Viking Road, Epping