


INVACARE MATRX LIBRA CUSHION



The next generation cushion - Superior comfort with a high level of protection and positioning


- Improved design combines foam and fluid to provide the highest level of skin protection, positioning and adjustability for unique user needs
- New dual layer fluid sac using invacare's **Flo-tech™** technology for enhanced shear management and pressure redistribution
- Specially contoured to provide balance and stability to the pelvis, optimizing body function
- Lightweight and maintenance free, the **Libra** balances both functional needs with an active lifestyle
- Recommended for all users including those with specific individual seating needs

Features and options




Immersion

Maximises surface contact area in order to reduce peak pressures. Pressure = Force/Area.



Off-Loading

Transfers forces away from the ischial tuberosities with greater weight bearing on the trochanters, hips and thighs.



Envelopment

Fluid sac completely surrounds the bony prominences, even during position changes.



Pelvic Stability

Contouring of the pelvic well provides lateral trochanteric support and anterior/posterior pelvic stability, promoting spinal extension for improved upper extremity function, head control and physiological function.

Technical data



Matrx Libra

Weight: 1,6 kg



Max User Weight:
Std 136 kg
HD 227kg



Outer cover:
machine wash
80°



Inner liner:
surface clean
with damp cloth

Matrx Libra Customisation Kit



Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za