

KHULA BUGGY

CODE:3410

Modern and sporty design
Full aluminium powder coated frame
Compact folding
Individual adjustability thanks to various adjustment options
Multifunctional seat unit
Configurable for unique individual care thanks to different accessories
Growing thanks to adjustable seat width and depth, back height and lower leg length
Transportable
Angle-adjustable and extension-type plastic footplate
8"*1.6" PU front all-direction casters
11"*2"PU rear all-direction wheels with wheel lock



SPECIFICATION

Name:	KHULA BUGGY
Load:	60kg
N.W:	37kg including standard accessories
Seat Width:	18-30cm
Seat Depth:	30-41cm
Recline:	90-157°
Tilt-in Space:	40°
Back Height:	58-76cm
Lower Leg Adj:	20-41cm
Castor:	8x1.6" solid castor
Rear Rim:	11x2" with PU tyre

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

KHULA BUGGY

CODE:3410



Accessories Available:

- Tray Table
- Canopy
- Various Chest
- Harness options
- Indoor Activity Frame



Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za