

KOGIA BATHLIFT

The *Aquatec Kogia* makes bathing simple and has been designed with **comfort, ease of use** and **safety** in mind. Weighing only 10kgs, the *Kogia* is **lightweight** yet able to support up to 140kgs. Having been tested against the toughest safety standards, the bathlifter has a robust and stable structure allowing you to bathe with peace of mind.

- Easy to use with a smooth riser mechanism
- Easy to assemble and maintain
- Lightweight and easy to move and handle



Technical data

	<i>Aquatec Kogia</i>
Seat width with side flaps folded in	370 mm
Seat width with side flaps folded out	715 mm
Seat depth	450 mm
Seat height at highest position	420 mm
Seat height at lowest position	60 mm
Backrest width	370 mm
Height of backrest	660 mm
Base area	570 x 266 mm
Total weight without battery	10 kg
Max. load	140 kg
Backrest can be reclined by up to	40°
Battery	Lead storage battery



Floating handset

This easy to use, floating hand control has been **ergonomically** designed and can be placed in its handset holder when not in use.



Open base plate

There is no need to remove the *Kogia* from the bath when you clean thanks to the open base plate.

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za