

MANUAL HI-LO HOSPITAL BED

Code: ALK06-A328P



4 sections
3 function attendant hand crank

Specifications:

Width:	98cm
Length:	213cm
Height:	48 – 72cm
Load Capacity:	150kg
Back Rise:	70 degrees
Foot Lift:	30 degrees at the knee
Material:	Powder Coated Mild Steel
Head & Foot:	Removable ABS Head & Foot Board
Castors:	5" steel hub rubber castor with brake
STD Accessories:	Cotsides, 12cm Taffeta covered Mattress

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za