

MATRIX MINI BACK

Transforming kids seating

The *Invacare Matrix Mini Back* is contoured to fit children of all shapes and sizes, combining exceptional support in an appealing low profile package. The *Mini Back* features Mini Set hardware, a specially designed miniature version of the popular Easy Set hardware. Every element of this back support has been carefully crafted and constructed, keeping pediatric users at the forefront at every stage.



Lightweight back shell

The durable aluminum shell has been designed for easy attachment of the optional headrest.

Water-resistant, breathable cover

Cover and foam overlap the shell for comfort and protection.

Lumbo Sacral Support

Technical data

Mini	 Height 18 cm, 25 cm, 30 cm & 36 cm	 Width 25 cm - 36 cm	 Depth* Adjustment 10 cm	 Angle Adjustment 40° range +/- 20°	 Product Weight 0,8 kg and up
Mini	 Weight Capacity 68 kg	 Limited Warranty 2 Years	* Limitation of depth/angle adjustment if used on chairs narrower than model width		

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za