

PORTABLE POOL HOIST

PRODUCT CODE: NART1000



The Nart Pool Lift is a State-of-the-art rechargeable battery powered pool lift system, stainless steel features and powder coating. Suitable for private pools as well as in a therapy use or recreational facility which has a pool . This pool lift provides users Easy pool access with an Easy set up as well as easy demount feature whenever you wish. The Nart Pool Lift can be used in outdoor and indoor pools. The Nart Pool Lift doesn't need any modification on the pools outside exterior. This lift can be easily removed without leaving any major hole or bulge on the floor. The Nart Pool Lift remains ready for use whenever you wish thanks to its external charging unit.

Technical Specifications

- Emergency button.
- Water resistant control.
- Detachable and rechargeable battery with external charging adapter
- Wheels suitable for wet surfaces
- Stroke length is 1250 mm. (610 mm above water 640 mm in water)
- Vacuum emergency stoppers
- Medical rustproof wheels with wheel safety chocks
- The device battery can be charged with 220 Volts.
- The device is 24 Volts
- Foldable armrest and footrest
- It has a lifting capacity of 135 kg.
- The weight of the device is 200 kg.
- It has a water-repellent seat suitable for wet use.
- Body is made of ST-52 steel and zinc plated and painted with electrostatic paint.
- The parts entering the water are 304 stainless steel and coated with electrostatic paint.
- Landing time to the pool is 45 seconds.
- The exit time from the pool is 55 second

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za