

MERITS TITAN POWER WHEELCHAIR

P200

Details:

Weight Capacity: 136 kg / 300 lb
Product Weight: 91 kg / 200 lb-95 kg / 209 lb
Driving Range: 40 km / 25 mi.
Max Speed: 7.2 km/h , 4.5 mph
Sizes: 16" & 18"

- Swing-in-out footrests provide no interference with surroundings.
- Adjustable seat angle 0°-16° / adjustable backrest angle 0°-30°.
- Height adjustable armrest and transparent side panel, comfortable cushion.
- Foldable backrest folds down for transport and storage.
- Wide armpad with central recess and bent end for better arm support.
- Rear suspension for better user experience.
- Removable battery box – Tool free removable battery box.
- Main chassis only 18kg (heaviest part).
- Optional lighting system and elevating legrests.



	16"
Driving Range	40 km / 25 mi.
Max Speed	7.2 km/h , 4.5 mph
Battery Type	12V / 22NF x 2 pcs
Motor	DC 24V, 180W,(Rating) / 630W(Max)
Max Climbing Angle	12°
Turning Radius	61 cm / 24"
Ground Clearance	10 cm / 4"
Controller	PG VS2 50
Charger	5A off-board
Weight Capacity	136 kg / 300 lb
Total Weight	91 kg / 200 lb
Front wheel	9" (2.8 / 2.5-4) foam filled tire
Drive Wheel	14" (3.00-8) foam filled tire
Rear wheel	14" (3.00-8) foam filled tire

	18"
Driving Range	40 km / 25 mi.
Max Speed	7.2 km/h , 4.5 mph
Battery Type	12V / 22NF x 2 pcs
Motor	DC 24V, 180W,(Rating) / 630W(Max)
Max Climbing Angle	12°
Turning Radius	61 cm / 24"
Ground Clearance	10 cm / 4"
Controller	PG VS2 50
Charger	5A off-board
Weight Capacity	136 kg / 300 lb
Total Weight	93 kg / 205 lb
Front wheel	9" (2.8 / 2.5-4) foam filled tire
Drive Wheel	14" (3.00-8) foam filled tire
Rear wheel	14" (3.00-8) foam filled tire