


SAGITTA SCOOTER

PRODUCT CODE: 300677
NAPPI CODE



Technical Specifications:

Motor	320 W
Controller	90 A
Seat depth	390 mm
Seat width	440 mm
Dimensions unfolded (L x W x H)	1255 x 560 x 965 mm
Max. weight of user	150kg
Max. speed (driving forward)	10 Km/h
Min. braking distance at max. speed	2100 mm
Angle seat surface	8°
Height of footplate	150 mm
Maximum height of obstacle	50 mm
Maximum safe slope	6°
Minimum turning radius	1500 mm
Minimum width required to reverse	1730 mm

	Version 1	Version 2
Battery	38 Ah	52 Ah
Battery charger	DC24V, 5Amp	DC24V, 6Amp
Total mass (with battery)	100 Kg	105 Kg
Total mass (w/o battery)	75 Kg	75 Kg
Maximum weight of single part	22 Kg	22 Kg
Endurance mileage	38 Km	50 Km