

# SPLIT II ECONOMY POWER

Y209



Frame:	Steel, powder coating
Controller:	Yanteon Controller
Motor:	250W*2, brush motor with electromagnetic brake
Battery:	Lead-acid 2 x 12V 20Ah
Charger:	HP 24V, 3A, Taiwan
Mas speed:	7 km/h
Travel Mileage:	15km
Front wheel:	8" PU caster
Rear wheel:	12.5" solid tire with steel spoke
Armrest:	Flip up slope armrest
Footrest:	Detachable, heel strap
Seat:	18"
Special:	Backrest half folding, Anti-tipper, Safety belt
Carton size:	87x48x81cm
G.W./N.W.:	37.5kg/35.5kg

Overall Length (cm)		108
Overall Width (cm)		63
Overall Height (cm)		97
Overall Width Folded		46
Seat Width (cm)		45
Seat Depth (cm)		44
Back Height (cm)		44
Floor to seat Height		53
Weight Capacity (kg)		135
Front Caster Spec.		8"
Rear Wheel Spec.		12.5"
G./N. Weight (kg)		37.0/35.0
Carton Size (cm)		87x38x81