

# STAR POWER WHEELCHAIR

PRODUCT CODE: STAR



The Star Electric Wheelchair is an impressive rear-wheel drive mobility solution, built with a durable powder-coated steel frame for strength and longevity. It comes equipped with a range of standard features designed for comfort and functionality, including power tilt and recline, a supportive headrest, integrated lights, height-adjustable armrests, and elevating leg rests. The seat and backrest are made from comfort injection foam, providing excellent support for extended use. Powered by two 340W motors and enhanced with rear suspension, the Star Electric Wheelchair delivers a smooth and powerful ride over various terrains.



## Product Specifications




Seat Width:  
43-53cm



Seat Depth:  
43-49cm



Seat Height:  
48cm



Front Wheel:  
9" PU Wheel



Rear Wheel:  
14" PU Wheel



L x W x H  
Dimensions:  
108x58x121cm



Gradeability:  
10 degrees



Driving Distance:  
30km\*



Max Speed:  
10km



Battery:  
24V 52AH



Unit Weight:  
106.5kg



Max Weight Capacity:  
135kg

[www.medop.co.za](http://www.medop.co.za) | [@medopcc](https://www.instagram.com/medopcc)



# Medop

Johannesburg:  
011 827 5893  
[enquiries@medop.co.za](mailto:enquiries@medop.co.za)  
Medopark, 1 Manchester Road,  
Wadeville

Cape Town:  
021 534 8628  
[admin@medop.co.za](mailto:admin@medop.co.za)  
Viking Business Park, Corner  
Park & Viking Road, Epping