

# STEED PRO TRANSFER CHAIR



CA6129-D

## ELECTRICAL TRANSFER CHAIR

- White powder coated frame
- Black backrest
- 5" front castors with brake & 3" rear wheels
- Electric lifting
- With bucket

CA6129

## MANUAL TRANSFER CHAIR

- White powder coated frame
- Black backrest
- 5" front castors with brake & 3" rear wheels
- Manual lifting
- With bucket

## MANUAL LIFTING PATIENT TRANSFER CHAIR

PRODUCT CODE : CA6129

- Feature:** Manual to lift or down  
**Material:** Steel metal frame, PE armrest, seat and backrest  
**Additional function:** brake on wheels, seat safety belt & buckle, feet support, commode seat design

### Specifications:

- Total height:** 83-102cm  
**Total length:** 78cm  
**Total width:** 60cm  
**Front Castor:** 5"  
**Rear Castor:** 3"  
**Load capacity:** 100kg  
**Packing size:** 70cmx29cmx70cm  
**Net weight** 28.4kg  
**Gross weight** 32kg

## ELECTRIC LIFTING PATIENT TRANSFER CHAIR

PRODUCT CODE: CA6129-D

- Feature:** Electrical remote control, press to lift or down  
**Material:** Steel metal frame, PE armrest, seat and backrest  
**Additional function:** brake on wheels, seat safety belt & buckle, feet support, commode seat design, hydraulic lifting labor-saving design, remote control

### Specifications:

- Total height:** 83-102cm  
**Total length:** 78cm  
**Total width:** 60cm  
**Front Castor:** 5"  
**Rear Castor:** 3"  
**Load capacity:** 100kg  
**Packing size:** 70cmx29cmx70cm  
**Net weight** 28.4kg  
**Gross weigh:** 32kg

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
 P.O. Box 51124, Raedene, 2124  
 enquiries@medop.co.za | www.medop.co.za  
 JHB Tel: 011 827 5893/4/5  
 CT Tel: 021 534 8628 | admin@medop.co.za