

Medop^{cc}

LIGHTWEIGHT FOLDING POWER WHEELCHAIR

CODE: DYW-459-46A8



- Aluminum powder coating frame, can be folded in small size;
- Flip-up armrest; foldable footplate;
- With pressure-reducing, washable, detachable "S" shape upholstery that are made by breathable material;
- Forward and backward adjustable controller can be installed in left armrest or right armrest freely;
- Motor: Brushless 180W×2 with Electromagnetic Brake;
- Range: 20±10% km;
- Max. speed: 6km/h;
- 8" PU casters; 8" solid rear wheel (Drive wheel);
- N.W.: 26.4KGs lithium battery included;
- Max. user weight: 100KGs ;
- Lithium battery: 24V/10Ah ;