

# WALKING TUTOR - HEAVY DUTY

- New aluminium oval tube frame construction.
- For institutions, rental or domestic use.
- Safe working load of 180kg, Australian standards certified.
- Fitted with new "light touch" hand brakes with a lock brake mechanism.
- The frame is widened between the rear wheels for maximum walking space.
- The Gutter arm walker is foldable for transportation.
- The oval tube provides a directional height adjustment function (Press button adjustment).
- The handles bars are multi adjustable for angle and fore/aft function.

**Registered design  
9210A**



W/Lock brake      Swivel castors

**9210A**

**9210B**

overall height	overall length	overall width	armrest length	folded width	castor size	unit weight	weight capacity
min103cm max127cm	88cm	65cm	min30cm max43cm	28cm	20cm	12.2kg	180kg

# Medop cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za