


VISION SPORT MID-WHEEL POWER WHEELCHAIR

P326A



- Mid-wheel drive design provides outstanding maneuverability in tight, compact spaces.
- Patented articulation mechanism allows front wheels climbing up the step/ threshold easily.
- Low center of gravity provides a smooth and stable ride both indoors and outdoors.
- Dual in-line motors for enhanced efficiency, torque, range and performance.
- Full suspension design offers a comfortable ride.
- Swing back controller bracket is standard.
- Height and width adjustable armrest, semi reclining backrest and height adjustable seat.

Driving Range	29 km / 18 mi.
Max Speed	8 km/h , 5 mph
Battery Type	12V U1x 2 pcs
Motor	DC 24V , 160W (Rating) / 520W (Max)
Max Climbing Angle	6°
Turning Radius	51 cm / 20"
Ground Clearance	5 cm / 2"
Controller	PG nVR2 50A / DY LiNX
Charger	4A off-board
Weight Capacity	136 kg / 300 lb
Total Weight	84 kg / 185 lb
Front wheel	F/R 6" x 2" PU tire
Drive Wheel	10" (260 x 85) Foam filled tire
Rear wheel	F/R 6" x 2" PU tire



Width 61cm	Seat Width 46cm	Seat Depth 41cm	Seat Height 55-60cm	Length 86cm	Footrest 43-53cm	Back Height 46cm	Height 123cm	Max Load 136kg	Product Weight 84kg
---------------	--------------------	--------------------	------------------------	----------------	---------------------	---------------------	-----------------	-------------------	------------------------

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za